

Colenta[®] ID-Printer TYPE B2N




Daylight ID-Printer System

(part# ID-B2N)

The Colenta ID Printer is for both AP and PA applications and incorporates the latest technology to record all the important patient data clearly and precisely onto each x-ray through the window area of the cassette in daylight operating conditions.




ID Card Dimensions



- Type.....AP/PA Daylight ID printing system used with window cassette
- Functions.....prints patient's data written or typed on the white card which user provides month, date and year with time are automatically printed
- Printing Method.....by 2 LED
- Required Cassette.....window type cassette
- Density Adjusting.....density is adjustable in 10 steps, indicated by display
- Safeguard.....The printer works when ID card is properly inserted into the card holder. 6 types error indication
- Clock.....built-in digital clock is battery backup. It prints month, date, year and time.
- Power Source.....single phase, AC 110-240V, 50/60Hz
- Measurements and Weight. .324 mm (L) x 435 mm (W) x 313 mm (H), about 6,5 kg

Printed sample on film



Printing dimensions are 22,5 mm x 73 mm - plus date and time.



Technical specification subject to change without notice.

COLENTA Labortechnik GmbH & Co. KG
 Neunkirchner Strasse 117
 A-2700 Wiener Neustadt
 Tel.: +43-2622-28311-0
 Fax: +43-2622-28311-7
 E-Mail: office@colenta.at
<http://www.colenta.at>

ACI-COLENTA INTERNATIONAL Div.
 Forge Close, Little End Rd,
 Eaton Socon, Cambs. UK
 Tel.: +44-1480-211118
 Fax: +44-1480-215236
 Fax: +44-870-1321352
 E-Mail: aci@patflaherty.com